GMD UNIVERSAL URINARY INCONTINENCE SLING-1012

The GMD Universal Sling™ should be placed using one or more of the GMD T-Slot reusable trocars

Product	Code
GMD Universal Urinary Incontinence Sling™-1012	1012
GMD T-Slot Reusable Trocars	

Trocar Type	Metal Handled T-Slot Trocars Product Codes	Thermoplastic Handled T-Slot Trocars Product Codes
Reusable Trocar Right Spiral	2041-T	2041-2T
Reusable Trocar Left Spiral	2042-T	2042-2T
Reusable Trocar Short Curve	2043-T	2043-2T
Reusable Trocar Long Curve	2044-T	2044-2T
Reusable Trocar Large Hook	2046-T	2046-2T
Reusable Trocar Right Spiral IO	2047-T	2047-2T
Reusable Trocar Left Spiral IO	2048-T	2048-2T
Reusable Trocar Short Curve XL	2049-T	2049-2T

INSTRUCTIONS FOR USE

Rx ONLY Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

ORIGINAL

DCP# P-10-100

Made in USA
 05-0088 Rev A

GMD Universal Sling™

Please read all information carefully: Failure to properly follow instructions may result in improper functioning of the device and lead to injury.

Important:

This package insert is designed to provide instructions for use for the Generic Medical Devices, Inc. (GMD) Universal Sling™-1012 . The GMD Universal Sling™-1012 should be placed using one or more of the GMD reusable trocars. The GMD Universal Sling™-1012 is supplied sterile and is intended for single patient use only. Carefully examine the package to verify that neither the package nor the sterilized contents have been damaged in shipment. DO NOT USE if the sterile barrier on product is damaged. Immediately return damaged product to GMD.

The design of the GMD Universal Sling™-1012 allows the physician to utilize devices for placement of a polypropylene mesh at the urethra or bladder neck.

THIS PRODUCT IS INTENDED FOR USE ONLY BY PHYSICIANS (1) WITH ADEQUATE TRAINING AND EXPERIENCE IN TREATMENT OF FEMALE STRESS URINARY INCONTINENCE AND (2) WHO ARE FAMILIAR WITH AND EXPERIENCED IN THE PLACEMENT OF NON-ABSORBABLE MESH SLINGS USING THE STANDARD OUTSIDE-IN AND INSIDEOUT SURGICAL TECHNIQUES AND APPROACHES. THE PHYSICIAN IS ADVISED TO CONSULT THE MEDICAL LITERATURE REGARDING SPECIFIC TECHNIQUES, COMPLICATIONS, AND HAZARDS ASSOCIATED WITH THE INTENDED PROCEDURES. THE PHYSICIAN SHOULD BE FAMILIAR WITH SURGICAL PROCEDURES.

These instructions are intended for general use of the device. Variations in use may occur in specific procedures due to individual technique and patient anatomy.

ORIGINAL

TABLE OF CONTENTS

1. Symbols and Definitions3	6. Warnings and Precautions8
2. Device Description3	7. Potential Complications8
3. Indications for Use5	8. How Supplied9
4. Instructions for Use5-7	9. Storage.....9
5. Contraindications7	10. Warranty9

1. Symbols and Definitions

The following symbols are used in conjunction with this device:

	Single use only, do not reuse
	Expiration date, use-before date
	Batch code
	Consult Instructions for Use
	This Product Contains No Detectable Latex
UPN	Product Number
	No DEHP


2. Device Description

The GMD Universal Sling™-1012 is a sterile, single use device for the treatment of female stress urinary incontinence. The GMD Universal Sling™-1012 is comprised of a polypropylene knitted mesh protected by a disposable low density polyethylene sheath with a disposable universal low density polyethylene bullet tip holding a PTFE coated, non-absorbable suture loop at each end for attachment of the sling to GMD's T-Slot Reusable Trocars. The method of placement and surgical approach chosen by the physician should be appropriate for the patient's diagnosis and anatomy. The physician should choose the appropriate trocar for the chosen sling placement method.

ORIGINAL


GMD Universal Sling™-1012

Detail of GMD Universal Sling™-1012

The GMD Universal Sling™ should be placed by using one or more of the GMD reusable trocars.

GMD Reusable Trocars:

Trocar Type	Metal Handled T-Slot Trocar Product Code	Thermoplastic Handled T-Slot Trocars Product Codes
A: Right Spiral	2041-T	2041-2T
B: Left Spiral	2042-T	2042-2T
C: Short Curve	2043-T	2043-2T
D: Long Curve	2044-T	2044-2T
E: Large Hook	2046-T	2046-2T
F: Right Spiral IO	2047-T	2047-2T
G: Left Spiral IO	2048-T	2048-2T
H: Short Curve XL	2049-T	2049-2T

3. Indications for Use

The GMD Universal Sling™-1012 is intended for use in women as a suburethral sling for the treatment of stress urinary incontinence (SUI) resulting from either urethral hypermobility and/or intrinsic sphincter deficiency.

The use of the device in males and children under 18 years of age is not supported by clinical studies.

4. Instructions for Use


A. Prior to Use

The bladder should be empty prior to initiating the use of this product. Ensure that the bladder, urethra, and other important landmarks are properly identified and that no pathology is present that would be contraindicated in using this product.

1. Choose the appropriate sling placement method and surgical based on the patient's diagnosis and anatomy.
2. Choose the GMD trocar(s) appropriate for the chosen surgical approach.
3. Prepare the skin of the abdomen, trunk, perineum and vaginal operative sites and drape the patient using standard surgical techniques.
4. Make appropriate incisions and develop a periurethral pocket at the midurethra and/or bladder neck to facilitate passage of the trocar without the need for using stylets or catheter guides.

B. Preparation and use of the GMD Universal Sling™-1012 – INSIDE-OUT APPROACHES –

1. Load the suture loop into the T-Slot of the chosen trocar and pull the loop in the T-Slot towards the trocar handle.


ORIGINAL

- Advance the trocar through the periurethral pocket with one hand, while tensioning the suture loop to hold it firmly in the T-Slot with the other hand until the T-Slot and suture loop protrude from the lower abdominal incision.

- Insert both ends of the sling assembly using the trocars as appropriate for the surgical placement method chosen such that the sling assembly when placed correctly forms a loop beneath the urethra.
- Remove the suture loops from the trocar and pull the suture loop out through the tissue until approximately 2" to 3" (5-7,5 cm) of the suture loop is outside the abdominal incision.
- Confirm the sheath and mesh are not twisted before tensioning the sling assembly. Adjust the sling assembly by pulling outwards on the sleeves, so that the sling material is centered below the urethra.


- Appropriately tension the sling assembly according to physician preference.
- Cut the sheath / sling ends on each side below the interface with the bullet tip to remove each bullet tip/suture loop subassembly and expose each sheath end.
- Grasp each of the exposed sheath ends independent of the sling with a clamp on each side and pull outwards with equal tension on each sheath for removal. Stabilizing the sheath and sling material under the urethra will facilitate the removal of the sheath during this removal phase.
- Verify the tension of the implanted sling and adjust, as necessary.
- Cut the distal ends of the sling so that the sling ends retract below the skin incisions.
- Close all incisions according to standardized methods.

C. Preparation and use of the GMD Universal Sling™ -1012 – OUTSIDE-IN APPROACHES –


- Advance the trocar through the lower abdominal incision until the T-Slot protrudes from the periurethral pocket.

2. While holding the trocar in the proper position with one hand place the suture loop into the T-Slot and pull the suture loop towards the trocar needle tip.

3. Pull the trocar and attached suture loop back through the tissue until approximately 2" to 3" (5-7,5 cm) of the suture loop is protruding from the lower abdominal incision. Remove the suture loop and sling assembly from the trocar.
4. Insert both ends of the sling assembly using the trocars as appropriate for the surgical placement method chosen such that the sling assembly when placed correctly forms a loop beneath the urethra.
5. Confirm the sheath and mesh are not twisted before tensioning the sling assembly. Adjust the sling assembly by pulling outwards on the suture loops, so that the sling material is centered below the urethra.

6. Appropriately tension the sling assembly according to physician preference.
7. Cut the sheath / sling ends on each side below the interface with the bullet tip to remove each suture loop/bullet tip subassembly and expose each sheath end.
8. Grasp each of the exposed sheath ends independent of the sling with a clamp on each side and pull outwards with equal tension on each sheath for removal. Stabilizing the sheath and sling material under the urethra will facilitate the removal of the sheath during this removal phase.
9. Verify the tension of the implanted sling and adjust, as necessary.
10. Cut the distal ends of the sling so that the sling ends retract below the skin incisions.
11. Close all incisions according to standardized methods.

5. Contraindications

The mesh suburethral sling implant is contraindicated in the following patients:

- This procedure should not be performed in pregnant patients.
- The procedure should not be performed in patients with future growth potential including women with future plans of pregnancy.
- This procedure should not be performed with patients who are on anti-coagulation therapy or patients who have a urinary tract infection.

ORIGINAL

6. Warnings and Precautions

- **The GMD Universal Sling™-1012 should be placed using one or more of the GMD single use or reusable trocars.**
- **Do not use** GMD Universal Sling™-1012 for patients who are on anti-coagulation therapy.
- **Do not use** GMD Universal Sling™-1012 for patients who have a urinary tract infection.
- Physicians using the GMD Universal Sling™-1012 should be familiar with and experienced in the use of mesh slings for the treatment of female urinary incontinence and the various surgical approaches for the placement of those slings.
- Acceptable surgical practice should be followed for the GMD Universal Sling™-1012 for the management of contaminated or infected wounds.
- The surgical procedure should be performed with care to avoid large vessels, nerves, bladder and bowel. Attention to patient anatomy and correct passage of the device will minimize risks.
- Bleeding may occur post-operatively. Observe for any symptoms or signs before releasing the patient from the surgical recovery process.
- Standard surgical practices should be followed for the sling procedure, as well as for the management of contaminated or infected wounds.
- The procedure should be performed with very careful attention to avoid laceration of any vessels, nerves, bladder, and bowel.
- Do not remove the protective sheath covering the mesh implant until proper position has been confirmed.
- Ensure the mesh is placed without tension under the urethra and/or bladder neck.
- Use of this device should be done with the understanding that subsequent infection may require removal of the mesh.
- Patients should be counseled to refrain from heavy lifting, exercise, and intercourse for a minimum of four (4) weeks after the procedure. The physician should determine when it is suitable for each patient to return to normal activities.
- If dysuria, bleeding, or other problems occur, the patient should be instructed to contact the physician immediately.
- Do not use any mechanical means of contact with the mesh (such as clips, staples, etc.) within the urethral support region of the mesh, as mechanical damage to the mesh may occur.
- Avoid excessive tension on the mesh during handling.

7. Potential Complications

The following complications have been reported due to sling placement, but are not limited to:

- As with all implants, local irritation at the wound site and/or a foreign body may occur.
- Tissue responses to the implant could include vaginal extrusion, erosion through the urethra or other surrounding tissue, migration of the device from the desired location, fistula formation, and inflammation. The occurrence of these responses may require removal of the entire mesh.
- Like all foreign bodies, the mesh may potentiate an existing infection.
- Excess tension may cause temporary or permanent lower urinary tract obstruction and retention.

- Known risks of surgical procedures for the treatment of incontinence include pain, infection, erosion, device migration, and complete failure of the procedure resulting in incontinence and mild to moderate incontinence due to incomplete support or overactive bladder.
- In addition to the above listed potential complications, allergic reaction, abscess, detrusor instability, pelvic and vaginal pain, dyspareunia, vaginal bleeding, vaginal discharge, dehiscence of vaginal incision, edema and erythema at the wound site have been reported due to a sling procedure.
- It has also been reported that groin pain, orthostatic symptoms, fatigue, and shortness of breath may occur due to potential development of hematoma in the pelvic area.

8. How Supplied

The GMD Universal Sling™-1012 is provided sterile (ethylene oxide) for single use. Do not resterilize. Do not use package if opened or damaged. Discard opened unused devices. If damage has occurred, contact your GMD representative.

After use, disposal of the product and packaging should be in accordance with hospital, administrative, and/or local government policy.

9. Storage

Store at controlled room temperature. Rotate inventory so that products are used prior the expiration date on package label. Do not use expired product.

10. Warranty

GMD warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, of the GMD Universal Sling™-1012, as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond GMD's control, directly affect the GMD Universal Sling™-1012 and the results obtained from its use. GMD's obligation under this warranty is limited to the repair or replacement of the product and GMD shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of the GMD Universal Sling™-1012. GMD neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **GMD assumes no liability with respect to instruments that are reused, reprocessed, or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments supplied in the GMD Universal Sling™.**

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All rights reserved
GMD, Inc.